


Introducing a new solution to an old problem

dispo'clean

unique

single passage cleaning process

for endoscopes



*Thoroughly & rapidly cleans all internal Endoscope channels
Suitable for use with all makes of Endoscopes*

*A revolution in traditional endoscope brushing techniques
Dispo'Clean offers, cost effective, safer, faster
and more efficacious cleaning than ever previously obtainable.*

EFFECTIVE

Dispo'Clean utilises an active group of three sets of flexible multi-blades to thoroughly clean all internal endoscope channels in one simple operation. The one piece can be used for all channels of the one scope after washing in an enzyme based cleaning solution between passes

EFFICIENT

A single one directional pass facilitates repeatable, faster and higher quality cleaning of lumens than traditional brushing

SAFETY

Dispo'Clean's unique proven patented process offers, safer, faster and more efficacious pre-disinfection cleaning

COMPREHENSIVE RANGE

Available in two sizes (colour coded) for use with all endoscopes with channel diameters between 2.8mm and 4.2mm

Choice of dispense formats:

- Box of 250 individually wrapped units
- Compact dispense reel of 250 units

CODE

DC2.8REEL	Green Gastro Cleaning Blades 2.8-3.2 Channel Size
DC2.8SINGLES	Green Gastro Cleaning Blades 2.8-3.2 Channel Size
DC3.2REEL	Yellow Colon Cleaning Blades 3.2-4.2 Channel Size
DC3.2SINGLES	Yellow Colon Cleaning Blades 3.2-4.2 Channel Size

DESCRIPTION

PACKAGING

Reel 250
Box 250 individually wrapped
Reel 250
Box 250 individually wrapped