

# FEBRUARY 2008 SPECIALS

Introducing

## WARMING CABINET/INCUBATOR

- Made entirely of inox AISI 304-18/10 stainless steel
- Heat insulated with fibreglass
- Electronic temperature control
- Temperature and cycle time display
- Temperature regulation between 30°C and 70°C
- Cycle time regulation from 0 to 120 min.
- Alarms for thermostat or resistances failures
- Temperature fall alarm
- Voltage 230V 50Hz option (110V or 60Hz)
- Door lockup with key
- Prearranged for wall fixing

Capacity: 104L  
Internal Measurements: 620 x 420 x 400 mm

Code	Description
CBM2429V	CBM WARMING CABINET 620X420X400

**PRICE: \$2,100.00**



## SILICONE MAT

- Sterilizable, fully washable and flexible
- Thermo-resistant (-40°C to 300°C)
- Suitable for central sterilization services
- Guaranteed protection for surgical instruments, trays, containers, micro-surgery and ultrasonic cleaners.

Code	Description
MSB500X320	SILICONE MAT BLUE 500 X 320 MM

**NORMAL WHOLESALE PRICE: \$240.93**

**SPECIAL PRICE: \$168.65**



## COLD INSTANT COMPRESS

- Easy activation
- Lasts around 15-40 minutes
- Cold compress can reach a low of -4°C
- Convenient and portable
- Size: 15 x 23 cm

Code	Description
INTICEUREA	INTCO INSTANT COLD COMPRESS (BOX 40)

**PRICE: \$76.25**



Special discount available for Fixed Term Contracts.

Please talk to your Account Manager.

**Cooper Medical**  
Quality Service

Special prices are valid until the end of February 2008 or while stock lasts. Prices exclude GST