

NEW PRODUCT INTRODUCTION

LIBERTY ULTRASONIC BATH

The affordable way to improve cleaning quality!

The unique and compact Liberty® Ultrasonic Bath provides effective and reliable cleaning in a convenient table-top design. For a surprisingly affordable price! It can thoroughly clean dirt and deposits from articles and effectively displaces contaminants. The Liberty® Ultrasonic Bath is ideal for surgeries, jewellers, dentists, laboratories and domestic use.

Liberty®

**NEW
PRODUCT**

Features:

- Stainless steel tank SUS304
- Includes basket and lid
- Tank Capacity: 1000 mL
- Power: 50 W
- Frequency: 40 KHz
- Weight: 1.6 kg

Measurements:

- Unit Size: 210 x 190 x 165 mm
- Tank Size ID: 165 x 135 x 70 mm
(Diagonal 190 mm)
- Basket Size ID: 160 x 125 x 56 mm
(Diagonal 185 mm)



Cooper Medical Price
\$260.00 ex GST

Code	Description
UBL	ULTRASONIC BATH

Cooper Medical
Quality Service