



CUSTOM MADE


Cooper Medical Ltd is a family owned and operated New Zealand business. We supply new and custom made medical equipment and surgical instruments. We source products from all over the world and from New Zealand manufacturers. We also provide Biomedical Engineering support.

Instrument stringers

Instrument stringers help keep instruments organised during sterilisation and this helps you maximize usable tray/container space. Made of high quality stainless steel.

Available in 3 sizes and additional sizes can be manufactured on request.

Code: STRINGER-10 Small: 100mm long x 63mm wide
Code: STRINGER-15 Medium: 150mm long x 63mm wide
Code: STRINGER-20 Large: 200mm long x 63mm wide



Sterilising boxes and trays with lids

Made with medical grade stainless steel 304. Sterilising box has closed ends. Perfect for smaller instruments.

Manufactured in New Zealand. Other sizes can be made on request.

Sizes:


- L 50mm x W 50mm x H 20mm
- L 100mm x W 50mm x H 20mm
- L 150mm x W 50mm x H 20mm
- L 210mm x W 50mm x H 20mm
- L 310mm x W 50mm x H 20mm



Stainless steel wire tray with lid

Large stainless steel wire tray with slot in lid. Customer used for bariatric retractors.

Length 750mm x width 240mm x height 60mm
Very durable, 304 grade medical stainless steel. Comes with a 5 year warranty.



www.coopermedical.co.nz

PH: 0800 0800 823 - Email: admin@coopermedical.co.nz - Sales - Murray Cooper Ph: 027 224 4374

57 Mahunga Drive, Mangere Bridge, Auckland 2022, New Zealand

Re-usable sterilisable quiver

Made of durable polyethelyne and teflon, with removable base. Manufactured in New Zealand.

Sizes:

- 200mm- order code:
- 250mm- order code:
- 310mm- order code:



Acrylic custom noise box / cleaning box

Made of durable clear acrylic, with removable base. Manufactured in New Zealand.

4.5mm clear acrylic- Front and back bent, as well as top bent in 50mm, with piano hinges for lid.

Front cut out is centered with an opening 300mm wide x 150 mm deep from the top edge.



Size: overall size 600mm wide x 400mm deep x 400mm high.

Lead apron hangers and trolley

Lead apron hangers, manufactured in New Zealand. Made of durable heavy duty stainless steel. Order code: LAH1001

R100 Mobile Lead apron trolley- fits x10 hangers. More information on request. Order code: LATR100



Upholstery

Upholstery in high quality vinyl for:

- Theatre squabs
- Lateral support pads
- Stools

