

THE ORIGINAL AND BEST



MAS-007 Single Use Cholangiogram Forceps

FEATURES

- Dimensions 5mm x 330mm
- Accepts up to 6fr Catheter

**Developed in Australia by
Michael Alan Swan (Endovision Founder)**

The MAS-007 continues to be the Gold Standard in Single Use Cholangiogram Forceps and is the preferred product by 90% of surgeons when compared to White Handled Single Use Cholangiogram Forceps.

Australian Designed. German Engineering Quality.

Endovision

ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT



Twelve consecutive Cholangiogram forceps, still functional, from five or six different manufacturers were dismantled and inspected for internal cleanliness. All twelve were contaminated as shown above. Since 2009, Endovision has offered the MAS-007 Single Use Cholangiogram Forceps.

Australian Designed. German Engineering Quality.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

03 9550 1819 | 0409 865 105 | Info@Endovision.com.au | Endovision.com.au